

FZD5 Rabbit pAb

货号: B12247

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:500
理论分子量	64kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat kidney
细胞定位	Cell membrane,Golgi apparatus membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-240 of human FZD5 (NP_003459.2).
序列	RMSCDRLPVLGRDAEVLCDYNRSEATTAPRPFPKPTLPGPPGAPASGGECPAGGPFVCKCREPFVPILKESHPLYNKV RTGQVPNCAVPCYQPSFSADERTFATFWIG

靶点信息

研究背景	Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD5 protein is believed to be the receptor for the Wnt5A ligand.
基因ID	7855

基因名	FZD5
Swiss	Q13467
别名	FZD5;C2orf31;HFZ5

产品验证



Western blot analysis of FZD5 expressed in Rat kidney using FZD5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn