

CEP164 Rabbit pAb

货号: B12260

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus IF: Mus musculus IHC: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	163kDa/164kDa
实测分子量	164kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375
细胞定位	Cytoplasm,Nucleus,centriole,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 730-829 of human CEP164 (NP_055771.4).
序列	EASEKSEQAALNAAKEKALQQLREQLEGERKEAVATLEKEHSAELERLCSSLEAKHREVVSSLQKKIQEAQQKEEAQLQKCLGQVEHRVHQSYHVAGYE

靶点信息

研究背景	This gene encodes a centrosomal protein involved in microtubule organization, DNA damage response, and chromosome segregation. The encoded protein is required for assembly of primary cilia and localizes to mature centrioles. Defects in this gene are a cause of nephronophthisis-related ciliopathies. Two transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	22897
基因名	CEP164
Swiss	Q9UPV0
别名	CEP164;NPHP15

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn