

USP26 Rabbit pAb

货号: B12287

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	104kDa
实测分子量	104kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U-937,Mouse testis,Mouse liver,Mouse heart,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 664-913 of human USP26 (NP_114113.1).
序列	TSLCQFHKAGGKPASSPGTPLSKVDFQTPENPKRKKYVKTSKFVAFDRIINPTKDLYEDKNIRIPERFQKVSEQTQQCDGMRICEQAPQQALPQSFPKPGTQGHTKNLLRPTKLNQKSNRNSLLALGSNKNPRNKDILDKIKSKAKETKRNDKGDHTYRLISVVSHLGKTLKSGHYICDAYDFEKQIWFTYDDMRVLGIQEAQMQEDRRCTGYIFFYMHNEIFEEMLKREENAQLNSKEVEETLQKE

靶点信息

研究背景	This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility.
基因ID	83844

基因名	USP26
Swiss	Q9BXU7
别名	USP26

产品验证



Western blot analysis of USP26 expressed in Mouse testis,Rat liver,HeLa using USP26 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn