

PIGR Rabbit pAb

货号: B12353

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus IF: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	83kDa
实测分子量	100KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Rat liver,Mouse lung,Mouse liver
细胞定位	Cell membrane,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 460-600 of human PIGR (NP_002635.2).
序列	GEPNLKVPGNVTAVLGETLKVPCHFCKFSSYEKYWCKWNNNTGCQALPSQDEGPSKAFVNCDENSRLVSLTLNLVTRAD EGWYWCGVKQGHFYGETAAVYVAVEERKAAGSRDVS LAKADAAPDEKVLDSGFREIENKAIQ

靶点信息

研究背景	This gene is a member of the immunoglobulin superfamily. The encoded poly-Ig receptor binds polymeric immunoglobulin molecules at the basolateral surface of epithelial cells; the complex is then transported across the cell to be secreted at the apical surface. A significant association was found between immunoglobulin A nephropathy and several SNPs in this gene.
------	--

基因ID	5284
基因名	PIGR
Swiss	P01833
别名	PIGR

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn