

LPTM5 Rabbit pAb

货号: B12366

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus Co-IP: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	
实测分子量	25KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,THP-1,U-937,Mouse spleen,Rat thymus
细胞定位	cytoplasmic vesicle,cytosol,lysosome,perinuclear region of cytoplasm,plasma membrane,transport vesicle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 153-262 of human LPTM5 (NP_006753.1).
序列	SSYMEVPTYLNFKSMNHMNYLPSQEDMPHNQFIKMMIIFSIAFITVLIFKVYMFKCVWRCYRLIKCMNSVEEKRNKMLQKV VLPSYEEALSLPSKTPEGGPAPPPYSEV

靶点信息

研究背景	This gene encodes a transmembrane receptor that is associated with lysosomes. The encoded protein, also known as E3 protein, may play a role in hematopoiesis. [provided by RefSeq, Feb 2009]
基因ID	7805

基因名	LAPTM5
Swiss	Q13571
别名	CLAST6;LAPTM5

产品验证



Western blot analysis of LAPTM5 expressed in Mouse spleen,Rat thymus,HL-60 using LAPTM5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn