

UBL3 Rabbit pAb

货号: B12405

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IHC: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200
理论分子量	13kDa
实测分子量	13kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,SW480,Mouse brain,Mouse kidney,Mouse lung,Rat brain
细胞定位	Cell membrane,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human UBL3 (NP_009037.1).
序列	MSSNVPADMINLRLLVSGKTKEFLFSPNDSASDIKHVYDNWPMDWEEEQVSSPNILRLIYQGRFLHGNVTLGALKLPFG KTTVMHLVARETLPEPNSQGQRNREKTGESNCCVIL

靶点信息

研究背景	Located in extracellular exosome.
基因ID	5412

基因名	UBL3
Swiss	O95164
别名	UBL3;HCG-1;PNSC1

产品验证



Western blot analysis of UBL3 expressed in Mouse brain,Rat brain,MCF7 using UBL3 Rabbit pAb at 1:100 0. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bloc king buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn