

MYH14 Rabbit pAb

货号: B12410

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Sus scrofa WB: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:1000 - 1:5000 IHC: 1:50 - 1:200
理论分子量	167kDa/202kDa/227kDa/228kDa/232kDa
实测分子量	252KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3
偶联物	Unconjugated
阳性对照	MCF7,U-937(Negative control),Rat lung
细胞定位	cytosol,extracellular exosome,growth cone
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide of human MYH14
序列	

靶点信息

研究背景	This gene encodes a member of the myosin superfamily. The protein represents a conventional non-muscle myosin; it should not be confused with the unconventional myosin-14 (MYO14). Myosins are actin-dependent motor proteins with diverse functions including regulation of cytokinesis, cell motility, and cell polarity. Mutations in this gene result in one form of autosomal dominant hearing impairment. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	79784

基因名	MYH14
Swiss	Q7Z406
别名	MYH14;DFNA4;DFNA4A;FP17425;MHC16;MYH17;NMHC II-C;NMHC-II-C;PNMHH;myosin;NMHCII-C;myosin-14

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn