

C1QTNF5 Rabbit pAb

货号: B12419

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	25kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,Mouse brain,Mouse kidney,Mouse uterus,Rat brain,Rat kidney
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 16-243 of human C1QTNF5 (NP_056460.1).
序列	SPPLDDNKIPSLCPGHPGLPGTPGHHGSQGLPGRDGRDGRDGAPGAPGEKGEGGRPGLPGRGDPGPRGEAGPAGPTGPAGECSVPPRSAFSAKRSESRVPPPSDAPLPFDRVLVNEQGHYDAVTGKFTCQVPGVYYFAVHATVYRASLQFDLVKNGESIASFFQFFGGWPKPASLSGGAMVRLEPEDQVWVQVGVDYIGIYASIKTDSTFSGFLVYSDWHSSPVFA

靶点信息

研究背景	This gene encodes a member of a family of proteins that function as components of basement membrane s and may play a role in cell adhesion. Mutations in this gene have been associated with late-onset retina l degeneration. The protein may be encoded by either a bicistronic transcript including sequence from th e upstream membrane frizzled-related protein gene (MFRP), or by a monocistronic transcript expressed fr om an internal promoter.
------	--

基因ID	114902
基因名	C1QTNF5
Swiss	Q9BXJ0
别名	C1QTNF5;CTRP5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn