

MEK6 Rabbit pAb

货号: B12433

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Oryza sativa , Sus scrofa
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	31kDa/37kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm,Nucleus,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human MEK6 (NP_002749.2).
序列	MSQSKGKKRNPGLKIPKEAFEQPQTSSTPPRDLDKACISIGNQNFEVKADDLEPIMELGRGAYGVVEKMRHVPSGQIMA VKRIRATVNSQEQKRLLMDLDISMRTVDCPFTVTFTYGALFREGDVWICMELMDTSLDKFYKQVIDKGQTIPEDILGKIAVSI VKALEHLHSLSVIHRDVKPSNVLINALGQVKMCDFGISGYLVDSVAKTIDAGCKPYMAPERINPELNQKGYSVKSDIWSL GITMIELAILRFPYDSWGTPFQQLKQVVEEPSQQLPADKFSAEFVDFTSQCLKKNSKERPTYPELMQHPFFTLHESKGTDVA SFVKLILGD

靶点信息

研究背景	This gene encodes a member of the dual specificity protein kinase family, which functions as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein phosphorylates and activates p38 MAP kinase in response to inflammatory cytokines or environmental stress. As an essential component of p38 MAP kinase mediated signal transduction pathway, this gene is involved in many cellular processes such as stress induced cell cycle arrest, transcription activation and apoptosis.
基因ID	5608
基因名	MAP2K6
Swiss	P52564
别名	MAP2K6;MAPKK6;MEK6;MKK6;PRKMK6;SAPKK-3;SAPKK3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn