

RAB11FIP5 Rabbit pAb

货号: B12444

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens Co-IP: Homo sapiens IF: Homo sapiens IP: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	70kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,PC-3
细胞定位	Cytoplasm,Cytoplasmic vesicle,Early endosome membrane,Golgi apparatus membrane,Mitochondrion membrane,Peripheral membrane protein,Recycling endosome membrane,secretory vesicle membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 354-653 of human RAB11FIP5 (NP_056285.1).
序列	HRSSISGSLPSSGSLQAVSSRFSEEGPRSTDDTWPRGSRSNSSEAVLGQEELSAQAKVLAPGASHPGEEEGARLPPEGKP VQVATPIVASSEAVAEGARKEERKPRMGLFHQQGLSRSELGRRSSLGEGGPILGASPHHSSSGEEKAKSSWFGLR EAKDPTQKPSHPVKPLSAAPVEGSPDRKQSRSSLSIALSSGLEKLTVTSGSIQPVTQAPQAGQMVDTKRLKDSAVLDQ SAKYHHLTHDELISLLLQRERELSQRDEHVQELESYIDRLLVRIMETSPTLLQIPPGPPK

靶点信息

研究背景	Rab effector involved in protein trafficking from apical recycling endosomes to the apical plasma membrane. Involved in insulin granule exocytosis. May regulate V-ATPase intracellular transport in response to extracellular acidosis.
基因ID	26056
基因名	RAB11FIP5
Swiss	Q9BXF6
别名	RAB11FIP5;GAF1;RIP11;pp75

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn