

HCK Rabbit pAb

货号: B12460

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IHC: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	57kDa/59kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60
细胞定位	Cell junction,Cell membrane,Cell projection,Cytoplasm,Cytoplasmic vesicle,Golgi apparatus,Lipid-anchor, Lysosome,Membrane,Nucleus,caveola,cytoskeleton,cytosol,cytosol,focal adhesion,podosome membrane, secretory vesicle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human HCK (N P_001165600.1).
序列	MGCMKSKFLQVGGNTFSKTETSASPHCPVYVPDPTSTIKPGPNSHNSNTPGIREAGSEDIIVVALYDYEAHHEDLSFQKG DQMVVLEESGEWWKARSLATRKEGYIPSNYVARVDSLETEEWWFFKGISRKDAERQLLAPGNMLGSFMIRDSETTKGSYSL SVRDYDPRQGDTVKH YKIRTLDNNGGFYISPRSTFSTLQELVDHYKKGNDGLCQKLSVPC

靶点信息

研究背景	The protein encoded by this gene is a member of the Src family of tyrosine kinases. This protein is primarily hemopoietic, particularly in cells of the myeloid and B-lymphoid lineages. It may help couple the Fc receptor to the activation of the respiratory burst. In addition, it may play a role in neutrophil migration and in the degranulation of neutrophils. Multiple isoforms with different subcellular distributions are produced due to both alternative splicing and the use of alternative translation initiation codons, including a non-AUG (CUG) codon.
基因ID	3055
基因名	HCK
Swiss	P08631
别名	HCK;JTK9;p59Hck;p61Hck

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn