

CD73/NT5E Rabbit pAb

货号: B12472

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus FC: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	57kDa/63kDa
实测分子量	70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,A375,Mouse kidney,Rat eye
细胞定位	Cell membrane,GPI-anchor,Lipid-anchor
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 32-264 of human CD73/NT5E (NP_002517.1).
序列	LHTNDVHSRLEQTSSESSKCVNASRCMGGVARLFTKVQQIRRAEPNVLLLDAGDQYQGTTWFTVYKGAEVAHFMNALRYDAMALGNHEFDNGVEGLIEPLLKEAKFPILSANIKAKGPLASQISGLYLPYKVLPGDEVVGIVGYTSKETPFLSNPGTNLVFDEITALQPEVDKCLKTLNVNKKIALGHSGFEMDKLIAQKVRGVDVVVGHSNTFLYTGNPPSKEVPAGKYP

靶点信息

研究背景	The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	4907
基因名	NT5E
Swiss	P21589
别名	NT5E;CALJA;CD73;E5NT;NT;NT5;NTE;eN;eNT

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn