

# C1QA Rabbit pAb

货号: B12485

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Mus musculus                                                                              |
| 应用    | <a href="#">WB</a>                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 26kDa                                                                                                |
| 实测分子量 | 25kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse lung,Mouse liver                                                                               |
| 细胞定位  | Secreted                                                                                             |
| 纯化    | Affinity purification                                                                                |


抗原信息

|      |                                                                                                                                                                                                                                       |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 23-245 of human C1QA (NP_057075.1).                                                                                                                     |
| 序列   | EDLCRAPDGKKGEAGRPGRGRPLKGEQGEPGAPGIRTGIQGLKGDQGEPSGNPGKVGYPGPGSGLGARGIPGIKG<br>TKGSPGNIKDQPRPAFSAIRRNPPMGGNVVIFDTVITNQEEPYQNHSGRFVCTVPGYYYFTFQVLSQWEICLSIVSSSRGQ<br>VRRSLGFCDTTNKGFLFQVVSAGGMVLQLQQGDQVWVEKDPKKGHIYQGSEADSVFSGFLIFPSA |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene encodes a major constituent of the human complement subcomponent C1q. C1q associates with C1r and C1s in order to yield the first component of the serum complement system. Deficiency of C1q has been associated with lupus erythematosus and glomerulonephritis. C1q is composed of 18 polypeptide chains: six A-chains, six B-chains, and six C-chains. Each chain contains a collagen-like region located near the N terminus and a C-terminal globular region. The A-, B-, and C-chains are arranged in the order A-C-B on chromosome 1. This gene encodes the A-chain polypeptide of human complement subcomponent C1q. |
| 基因ID  | 712                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 基因名   | C1QA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Swiss | P02745                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 别名    | C1QA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)