

VIP Rabbit pAb

货号: B12490

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus IHC: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,U-87MG,SH-SY5Y,Mouse spleen,Mouse brain,Rat brain
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-170 of human VIP (NP_003372.1).
序列	APSALRLGDRIPFEGANEPDQVSLKEDIDMLQNALAENDTPYYDVSRNARHADGVFTSDFSLLGQLSAKKYLESLMGKR VSSNISEDPVPVKRHSDAVFTDNYTRLRKQMAVKKYLNLSILNGKRSSEGESPDFPEELEK

靶点信息

研究背景	The protein encoded by this gene belongs to the glucagon family. It stimulates myocardial contractility, causes vasodilation, increases glycogenolysis, lowers arterial blood pressure and relaxes the smooth muscle of trachea, stomach and gall bladder. The protein also acts as an antimicrobial peptide with antibacterial and antifungal activity. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.
------	---

基因ID	7432
基因名	VIP
Swiss	P01282
别名	VIP;PHM27

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn