

MX1 Rabbit pAb

货号: B12492

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Dragon Li cats , Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	55kDa/75kDa
实测分子量	76KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Daudi
细胞定位	Cytoplasm,Cytoplasmic side,Endoplasmic reticulum membrane,Nucleus,Peripheral membrane protein,per inuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 368-662 of human MX1 (NP_001171517.1).
序列	ENEKMFFLIDKVNAFNQDITALMQGEETVGEEDIRLFTRLRHEFHKWSTIIENNFEQEGHKILSRKIQKFENQYRGRELPGFVN YRTFETIVKQIKALEEPAVDMLHTVTDMVRLAFTDVSIKNFEEFFNLHRTAKSKIEDIRAEQEREKEKLIRLHFQMEQIVYCQ DQVYRGALQKVREKELEEEKKKKSWDFGAFQSSSATDSSMEEIFQHLMAYHQEASKRISSHIPIIQFFMLQTYGQQLQKA MLQLLQDKDTYSWLLKERSDTSKDRKFLKERLARLTQARRRLAQFPG

靶点信息

研究背景	This gene encodes a guanosine triphosphate (GTP)-metabolizing protein that participates in the cellular antiviral response. The encoded protein is induced by type I and type II interferons and antagonizes the replication process of several different RNA and DNA viruses. There is a related gene located adjacent to this gene on chromosome 21, and there are multiple pseudogenes located in a cluster on chromosome 4. Alternative splicing results in multiple transcript variants.
基因ID	4599
基因名	MX1
Swiss	P20591
别名	MX1;IFI-78K;IFI78;MX;MxA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn