

SUN1 Rabbit pAb

货号: B12509

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	90kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	cytoplasm,nuclear envelope,nuclear membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SUN1 (O94901).
序列	MDFSRLHMYSPQCVPEPNTGYTYALSSSYSSDALDFETEHKLDPVFDSRMSRRSLRLATTACTLGDGEAVGADSGTSS AVSLKNRAARTTKQRRSTNKSAFSINHVSQVTSSGVSHGGTVSLQDAVTRRPPVLDESWIREQTTVDHFWGLDDDDGD LKGGNKAAIQNGDVGAAAATAHNGFSCSNCSMLSERKDVLTAHPAAPGPVSRVYSRDRNQKC

靶点信息

研究背景	This gene is a member of the unc-84 homolog family and encodes a nuclear nuclear envelope protein with an Unc84 (SUN) domain. The protein is involved in nuclear anchorage and migration. Alternatively spliced transcript variants have been described.
基因ID	23353

基因名	SUN1
Swiss	O94901
别名	SUN1;UNC84A

产品验证



Western blot analysis of SUN1 expressed in Mouse testis using SUN1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn