

DCAF11 Rabbit pAb

货号: B12519

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB
推荐浓度	WB: 1:1000 - 1:5000
理论分子量	50kDa/58kDa/61kDa
实测分子量	61KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human DCAF1 1 (NP_079506.3).
序列	MGSRNSSSAGSGSGDPSEGLPRRGAGLRRSEEEEEDEDVDLAQVLAYLLRRGQVRLVQGGGAANLQFIQALLDSEEN DRAWDGRLGDRYNPPVDATPDTRLEFNEIKTQVELATGQLGLRRAAQKHSFPRMLHQRERGLCHRGSFSLGEQSRVIS HFLPNDLGFTDS

靶点信息

研究背景	This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates. Multiple alternatively spliced transcript variants have been found for this gene.
基因ID	80344

基因名	DCAF11
Swiss	Q8TEB1
别名	DCAF11;GL014;PRO2389;WDR23

产品验证



Western blot analysis of DCAF11 expressed in Mouse liver using DCAF11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn