

MFN1 Rabbit pAb

货号: B12520

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Sus scrofa , Gallus gallus
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	83kDa
实测分子量	84kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse kidney,Rat kidney
细胞定位	Cytoplasm,Mitochondrion outer membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MFN1 (NP_077162.2).
序列	MAETVSPLKHFVLAKKAITAIFGQLLEFVTEGSHFVEATYRNPELDRIASEDDLVEIQGYRNKLAVIGEVLSRRHMKVAFFGR TSSGKSSVINAMLWDKV

靶点信息

研究背景	The protein encoded by this gene is a mediator of mitochondrial fusion. This protein and mitofusin 2 are homologs of the Drosophila protein fuzzy onion (Fzo). They are mitochondrial membrane proteins that interact with each other to facilitate mitochondrial targeting.
基因ID	67414

基因名	
Swiss	
别名	MFN1;hfzo1;hfzo2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn