

RSPH14 Rabbit pAb

货号: B12521

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	38kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,LO2,293T,Raji,Mouse brain,Rat brain
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-348 of human RSPH14 (NP_055248.1).
序列	MAHSQNSLELPININATQITTAYGHRALPKLKEELQSEDLQTRQKALMALCDLMHDPECIYKAMNIGCMENLKALLKDSNS MVRIKTTEVLHITASHVGRYAFLEHDIVLALSFLNDPSPVCRGNLYKAYMQLVQVPRGAQEIIISKGLISSLVWKLQVEVEE EEFQEFILDTLVLCLQEDATEALGSNVVLVLKQLLSANQNIRSKAARALLNVSISREGKKQVCHFVIPILVHLLKDPVEHVK SNAAGALMFATVITEGKYAALEAQAIGLLLELLHSPMTIARLNATKALTMLAEAEGRKALQTHVPTFRAMEVETYEKPQVAE ALQRAARIAISVIEFKP

靶点信息

研究背景	This gene encodes a protein with no known function but with slight similarity to a yeast vacuolar protein. The gene is located in a region deleted in pediatric rhabdoid tumors of the brain, kidney and soft tissues, but mutations in this gene have not been associated with the disease.
------	--

基因ID	27156
基因名	RSPH14
Swiss	Q9UHP6
别名	RSPH14;RTDR1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn