

TRAM1 Rabbit pAb

货号: B12529

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Macaca mulatta , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa/43kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,mouse liver
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 305-374 of human TRA M1 (NP_055109.1).
序列	CVTQAFMMWKFINFQLRRWREHSAFQAPAVKKKPTVTKGRSSKKGTEGVNGTLTSNVADSPRNKKEKSS

靶点信息

研究背景	This gene encodes a multi-pass membrane protein that is part of the mammalian endoplasmic reticulum. The encoded protein influences glycosylation and facilitates the translocation of secretory proteins across the endoplasmic reticulum membrane by regulating which domains of the nascent polypeptide chain are visible to the cytosol during a translocational pause.
基因ID	23471

基因名	TRAM1
Swiss	Q15629
别名	TRAM1;PNAS8;TRAM;TRAMP

产品验证



Western blot analysis of TRAM1 expressed in U-251MG,mouse liver using TRAM1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn