

KNTC1 Rabbit pAb

货号: **B12584**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa/250kDa
实测分子量	250kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Chromosome,Cytoplasm,Nucleus,centromere,cytoskeleton,kinetochore,spindle
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 2084-2209 of human KNTC1 (NP_055523.1).
序列	HSEKRHHQIKNFLGSCDPQVILKQLEEHMNTGQLAGFSHQIRSLLNININKKEFGILAKTKYFQMLKMHAMNTNNITELVNYLANDLSLDEASVLITEYSKHCGKPVPPDTAPCEILKMFLSGLS

靶点信息

研究背景	This gene encodes a protein that is one of many involved in mechanisms to ensure proper chromosome segregation during cell division. Experimental evidence indicated that the encoded protein functioned in a similar manner to that of the Drosophila rough deal protein.
基因ID	9735

基因名	KNTC1
Swiss	P50748
别名	KNTC1;ROD

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn