

# DECR1 Rabbit pAb

货号: **B12587**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Rattus norvegicus , Homo sapiens <b>IHC:</b> Rattus norvegicus
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	34kDa/36kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,BxPC-3,293T,A-549,Mouse brain,Mouse liver,Rat liver
细胞定位	Mitochondrion
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 35-335 of human DECR 1 (NP_001350.1).
序列	NTEALQSKFFSPLQKAMLPPNSFQGKVAFITGGGTGLGKGMTLLSSLGAQCVIASRKMDVLKATAEQISSQTGNKVHAI QCDVRDPDMVQNTVSELIKVAGHPNIVINNAAGNFISPTERLSPNAWKTTIDIVLNGTAFVTLLEIGKQLIKAQKGAFLSITTI YAETGSGFVPSASAKAGVEAMSKSLAAEWGKYGMRFNVIQPGPIKTKGAFSRLDPTGTTFEKEMIGRIPCGRGLTVEELAN LAAFLCSDYASWINGAVIKFDGGEEVLISGEFNDLRKVTKEQWDTIEELIRKTKGS

## 靶点信息

研究背景	This gene encodes an accessory enzyme which participates in the beta-oxidation and metabolism of unsaturated fatty enoyl-CoA esters.
基因ID	1666

基因名	DEC1
Swiss	Q16698
别名	DEC1;DEC1;NADPH;SDR18C1;2

## 产品验证



Western blot analysis of DEC1 expressed in Mouse brain,Rat liver,HepG2 using DEC1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)