

MFGE8 Rabbit pAb

货号: **B12612**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens,Mus musculus IF: Homo sapiens,Mus musculus WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	35kDa/37kDa/43kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	Membrane,Peripheral membrane protein,Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human MFGE8 (NP_001108086.1).
序列	MPRPRLLAALCGALLCAPSLLVALDICKNPCHNGGLCEEISQEVRGDVFPSTCTCLKGYAGNHCETKCVEPLGLENGNIA NSQIAASSVRVTFLGLQHWVPELARLNLRAGMVNAWTPSSNDDNPWIQVNLRRMWVTGVVTQGASRLASHEYLKAFKV AYSLNGHEFDIFIHDVNKKHKEFVGWNWKNNAVHVNLFETPV

靶点信息

研究背景	This gene encodes a preproprotein that is proteolytically processed to form multiple protein products. The major encoded protein product, lactadherin, is a membrane glycoprotein that promotes phagocytosis of apoptotic cells. This protein has also been implicated in wound healing, autoimmune disease, and cancer. Lactadherin can be further processed to form a smaller cleavage product, medin, which comprises the major protein component of aortic medial amyloid (AMA). Alternative splicing results in multiple transcript variants.
基因ID	4240
基因名	MFGE8
Swiss	Q08431
别名	MFGE8;BA46;EDIL1;HMFG;HsT19888;MFG-E8;MFGM;OAcGD3S;SED1;SPAG10;hP47

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn