

MICU2 Rabbit pAb

货号: B12616

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus , Homo sapiens IP: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	49kDa
实测分子量	45KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-87MG
细胞定位	Mitochondrion intermembrane space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human MICU2 (NP_689939.1).
序列	AWHHSRVSAARDGSFTVSAQKNVEHGIIYGKPSLRKQRFMQFSSLEHEGEYYMTPRDFLFSVMFEQMERKTSVKKLTK KDIEDTLSGIQTAGCGSTFFR

靶点信息

研究背景	Key regulator of mitochondrial calcium uniporter (MCU required to limit calcium uptake by MCU when cytoplasmic calcium is low. MICU1 and MICU2 form a disulfide-linked heterodimer that stimulate and inhibit MCU activity, depending on the concentration of calcium. MICU2 acts as a gatekeeper of MCU that senses calcium level via its EF-hand domains: prevents channel opening at resting calcium, avoiding energy dissipation and cell-death triggering.
------	---

基因ID	221154
基因名	MICU2
Swiss	Q8IYU8
别名	MICU2;1110008L20Rik;EFHA1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn