

ZIP14 Rabbit pAb

货号: B12621

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Pelteobagrus fulvidraco , Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	52kDa/54kDa
实测分子量	54KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	Cell membrane,Cell projection,Multi-pass membrane protein,lamellipodium
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-157 of human ZIP14 (NP_001121903.1).
序列	SSLGAPAIASAASFLQDLIHRYGEGDSLTLQQLKALLNHLDVGVGRGNVTQHVQGHRNLSTCFSSGDLFTAHNFSEQSRIG SSELQEFCPTILQQLDSRACTSENQENEENEQTEEGRPSAVEVWGYG

靶点信息

研究背景	Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A14 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008]
------	--

基因ID	23516
基因名	SLC39A14
Swiss	Q15043
别名	SLC39A14;HMNDYT2;LZT-Hs4;NET34;ZIP14;cig19

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn