

SLC6A14 Rabbit pAb

货号: B12651

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	72kDa
实测分子量	72KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,MCF7,BxPC-3,Mouse large intestine,Mouse lung,Mouse testis
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 131-234 of human SLC6 A14 (NP_009162.1).
序列	YYNVIIAYSLYMFASFQSELPWKNCSWSDKNCSRSPIVTHCNVSTVNKGIQEIIMNKSWVDINNFTCINGSEIYQPGQL PSEQYWNKVALQRSSGMNETGV

靶点信息

研究背景	This gene encodes a member of the solute carrier family 6. Members of this family are sodium and chlori de dependent neurotransmitter transporters. The encoded protein transports both neutral and cationic a mino acids. This protein may also function as a beta-alanine carrier. Mutations in this gene may be associ ated with X-linked obesity. A pseudogene of this gene is found on chromosome X.
基因ID	11254

基因名	SLC6A14
Swiss	Q9UN76
别名	SLC6A14;BMIQ11

产品验证



Western blot analysis of SLC6A14 expressed in Mouse large intestine,HeLa,HepG2 using SLC6A14 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn