

PIK3CG Rabbit pAb

货号: B12694

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	126kDa
实测分子量	126kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cell membrane,Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human PIK3CG (NP_002640.2).
序列	MELENYKQPVVLREDNCRRRRRMKPRSAASLSSMELIPIEFVLPTSQRKCKSPETALLHVAGHGNVEQMKAQVWLRALET SVAADFYHRLGPHHFLLLYQKKGQWYEIYDKYQVVQTLDCRLYWKATHRSPGQIHLVQRHPPSEESQAFQRQLTALIGYD VTDVSNVHDDLEFTRRGLVTPRMAEVASRDPKLYAMHPWVTSKPLPEYLWKKIANNCIFIVIHRSTTSQTIKVSPDDTPGA ILQSFFTKMAKKKSLMDIPESQSEQDFVLRVCGRDEYLVGETPIKNFQWVRHCLKNGEEIHVVLDTP

靶点信息

研究背景	Phosphoinositide 3-kinases (PI3Ks) phosphorylate inositol lipids and are involved in the immune response. The protein encoded by this gene is a class I catalytic subunit of PI3K. Like other class I catalytic subunits (p110-alpha p110-beta, and p110-delta), the encoded protein binds a p85 regulatory subunit to form PI3K . This gene is located in a commonly deleted segment of chromosome 7 previously identified in myeloid le ukemias. Several transcript variants encoding the same protein have been found for this gene.
------	---

基因ID	5294
基因名	PIK3CG
Swiss	P48736
别名	PIK3CG;PI3CG;PI3K;PI3Kgamma;PIK3;p110gamma;p120-PI3K

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn