

GALP Rabbit pAb

货号: B12727

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Sus scrofa , Rattus norvegicus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	13kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG
细胞定位	extracellular region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-116 of human GALP (NP_149097.1).
序列	APAHRGRGGWTLNSAGYLLGPVLHLPQMGDQDGKRETALEILDWLKAIDGLPYSHPPQPSKRNVMETFAKPEIGDLGML SMKIPKEEDVLKS

靶点信息

研究背景	This gene encodes a member of the galanin family of neuropeptides. The encoded protein binds galanin receptors 1, 2 and 3 with the highest affinity for galanin receptor 3 and has been implicated in biological processes involving the central nervous system including hypothalamic regulation of metabolism and reproduction. A peptide encoded by a splice variant of this gene, termed alarin, has vasoactive properties, displays antimicrobial activity against E. coli, and may serve as a marker for neuroblastic tumors.
基因ID	85569

基因名	GALP
Swiss	Q9UBC7
别名	GALP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn