

# MT-ATP6 Rabbit pAb

货号: B12731

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Homo sapiens , Mus musculus, Homo sapiens                                                 |
| 应用    | <a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a>                                        |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:100<br><b>IF/ICC:</b> 1:50 - 1:100                 |
| 理论分子量 | 25kDa                                                                                                |
| 实测分子量 | 25kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa,U-251MG,mouse brain,mouse heart,mouse kidney,rat kidney                                         |
| 细胞定位  | mitochondrial inner membrane,mitochondrial proton-transporting ATP synthase complex                  |
| 纯化    | Affinity purification                                                                                |


抗原信息

|      |                                                                                                             |
|------|-------------------------------------------------------------------------------------------------------------|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 1-50 of human MT-ATP6 (YP_003024 031.1). |
| 序列   | MNENLFASFIAPTILGLPAAVLILFPPLLIPTSKYLINNRLITTQWLI                                                            |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | Mitochondrial membrane ATP synthase (F1F0 ATP synthase or Complex V produces ATP from ADP in the presence of a proton gradient across the membrane which is generated by electron transport complexes of the respiratory chain. F-type ATPases consist of two structural domains, F(1 - containing the extramembraneous catalytic core and F(0 - containing the membrane proton channel, linked together by a central stalk and a peripheral stalk. During catalysis, ATP synthesis in the catalytic domain of F(1 is coupled via a rotary mechanism of the central stalk subunits to proton translocation. Key component of the proton channel; it may play a direct role in the translocation of protons across the membrane. |
| 基因ID  | 4508                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 基因名   | MT-ATP6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Swiss | P00846                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 别名    | ATPase6;MTATP6;MT-ATP6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)