

OsRDR6 Rabbit pAb

货号: B12855

产品信息

反应	Oryza sativa
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Oryza sativa
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	137kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	stems (before flower),leaves sheath (before flower),flower buds (before flower),flag leaves (after flower),flowers,milk seeds
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of Oryza sativa OsRDR6.
序列	

靶点信息

研究背景	nvolved in the RNA silencing pathway. Probably required for the generation of small interfering RNAs (siRNAs). Regulates shoot apical meristem (SAM) initiation and maintenance and leaf polarization through the trans-acting siRNAs (ta-siRNAs) pathway which probably modulates the expression of the ARF2, ARF3, ARF4, ARF14 and ARF15 genes.
基因ID	4323938

基因名	
Swiss	Q8LHH9
别名	SHL2; SHOOTLESS2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn