

PTPIP51/RMDN3 Rabbit pAb

货号: B12908

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mouse heart,Rat heart , Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	38kDa/52kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,HeLa,MCF7,Mouse liver,Mouse brain,Mouse kidney
细胞定位	Cytoplasm,Mitochondrion membrane,Mitochondrion outer membrane,Nucleus,Single-pass membrane protein,cytoskeleton,spindle,spindle pole
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 171-470 of human PTPIP51/RMDN3 (NP_060615.1).
序列	ESEGGYTTANAESDNERDSDKESEDGEDEVSCETVKMGRKDSLDEEEAASGASSALEAGSSGLEDVLP LLQQADELH RGDEQ GKREGFQLLLNNKLVYGS RQDFLWRLARAYS DMCELTEEVSEKSYALDGKEEAEEALEKGD ESADCHLWYAVL CGQLAEHESIQRRIQSGFSFKEHVDKAIALQPENPM AHFLLGRWCYQVSHLSWLEKKTATALLESPLSATVEDALQSFLKA EELQPGFSKAGRVYISKCYRELGKNSEARWWMKLALELPDVT KEDLAIQKDLEELEVILRD

靶点信息

研究背景	RMDN3 (regulator of microtubule dynamics 3) is a protein-coding gene. Diseases associated with RMDN3 include cerebritis, and noma. GO annotations related to this gene include protein binding. An important p aralog of this gene is RMDN1.
基因ID	55177
基因名	RMDN3
Swiss	Q96TC7
别名	RMDN3;FAM82A2;FAM82C;RMD-3;RMD3;ptpip51

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn