

# HACE1 Rabbit pAb

货号: B12918

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>IF:</b> Mouse bv2<br><b>WB:</b> Mus musculus,Homo sapiens                                         |
| 应用    | WB IHC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200                                                |
| 理论分子量 | 35kDa/63kDa/77kDa/102kDa                                                                             |
| 实测分子量 | 102kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | 293T,HeLa,Mouse brain,Mouse lung,Mouse spleen,Mouse testis                                           |
| 细胞定位  | Cytoplasm,Endoplasmic reticulum,Golgi apparatus,Golgi stack membrane                                 |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HACE1 (NP_065822.2). |
| 序列   | MERAMEQLNRLTRSLRRARTVELPEDNETAVYTLMPVMADQHRVSSELLSNSKFDVNYAFGRVKRSLHIAANCGSVECLVLLKKGANPNYQDISGCT      |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                             |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a HECT domain and ankyrin repeat-containing ubiquitin ligase. The encoded protein is involved in specific tagging of target proteins, leading to their subcellular localization or proteasomal degradation. The protein is a potential tumor suppressor and is involved in the pathophysiology of several tumors, including Wilm's tumor. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |             |
|-------|-------------|
| 基因ID  | 57531       |
| 基因名   | HACE1       |
| Swiss | Q8IYU2      |
| 别名    | HACE1;SPPRS |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)