

IL33 Rabbit pAb

货号: B12923

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Human corneal epithelial cells , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	16kDa/25kDa/26kDa/30kDa
实测分子量	33KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,A-549,Rat lung,Mouse lung,Mouse kidney,Mouse stomach
细胞定位	Chromosome,Cytoplasmic vesicle,Nucleus,Secreted,secretory vesicle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human IL33 (NP_254274.1).
序列	DSSITGISPITEYLASLSTYNDQSITFALEDESYEIYVEDLKKDEKKDKVLLSYYESQHPSPNESGDGVDGKMLMVTLSPTK

靶点信息

研究背景	The protein encoded by this gene is a cytokine that binds to the IL1RL1/ST2 receptor. The encoded protein is involved in the maturation of Th2 cells and the activation of mast cells, basophils, eosinophils and natural killer cells. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	90865
基因名	IL33

Swiss	O95760
别名	IL33;C9orf26;DVS27;IL1F11;NF-HEV;NFEHEV

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn