

KDM5B Rabbit pAb

货号: **B12924**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Human hepatocellular carcinoma cell line , Homo sapiens IP: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	175kDa/179kDa
实测分子量	176KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 784-883 of human KDM 5B (NP_006609.3).
序列	EESEMKKFPDNDLLRHLRLVTQDAEKCASVAQQLNGKRQTRYRSGGGKSQNQLTVNELRQFVTQLYALPCVLSQTPLLK DLLNRVEDFQHSQKLLSEE

靶点信息

研究背景	This gene encodes a lysine-specific histone demethylase that belongs to the jumonji/ARID domain-containing family of histone demethylases. The encoded protein is capable of demethylating tri-, di- and mono methylated lysine 4 of histone H3. This protein plays a role in the transcriptional repression or certain tum or suppressor genes and is upregulated in certain cancer cells. This protein may also play a role in genom e stability and DNA repair. Alternate splicing results in multiple transcript variants.
------	--

基因ID	10765
基因名	KDM5B
Swiss	Q9UGL1
别名	KDM5B;CT31;JARID1B;PLU-1;PLU1;PPP1R98;PUT1;RBBP2H1A;RBP2-H1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn