

U2AF1 Rabbit pAb

货号: B12964

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: 293T , Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	8kDa/19kDa/27kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,Mouse spleen
细胞定位	Nucleus,Nucleus speckle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human U2AF1 (NP_006749.1).
序列	MAEYLASIFGTEKDKVNCSFYFKIGACRHGDRCSRLHNKPTFSQTIALLNRYRNPQNSSQSADGLRCAVSDVEMQEHYDEF FEEVFTEMEKEYGEVEEMNVCDNLGDHLVGNVYVKFRREEDAEEKAVIDLNNRWFNQGPIHAELSPVTDFREACCRQYEM GECTRGGFCNFMHLKPISRELRRRELYGRRRKKHRSRSRSRERRSRSRDRGRGGGGGGGGGGGGGRERDRRRSRDRERS GRF

靶点信息

研究背景	This gene belongs to the splicing factor SR family of genes. U2 auxiliary factor, comprising a large and a small subunit, is a non-snRNP protein required for the binding of U2 snRNP to the pre-mRNA branch site. This gene encodes the small subunit which plays a critical role in both constitutive and enhancer-dependent RNA splicing by directly mediating interactions between the large subunit and proteins bound to the enhancers. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	7307
基因名	
Swiss	Q01081
别名	U2AF1;FP793;RN;RNU2AF1;U2AF35;U2AFBP

产品验证



Western blot analysis of U2AF1 expressed in 22Rv1, Mouse spleen using U2AF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn