

Progesterone Receptor Rabbit pAb

货号: B12968

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Ishikawa , Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	36kDa/38kDa/82kDa/87kDa/98kDa
实测分子量	118KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Mouse brain,Rat brain
细胞定位	Cytoplasm,Cytoplasm,Mitochondrion outer membrane,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Progest erone Receptor (NP_000917.3).
序列	MTELKAKGPRAPHVAGGPPSPEVGSPLLCRPAAGPFGSQTSDLPEVSAIPISLDGLLFPRPCQGQDPSDEKTQDQQSL SDVEGAYSRAEATRGAGGSSSSPPEKDSGLLDVLDTLAPSGPGQSQPSPACEVTSSWCLFGPELPEDPPAAPATQRV LSPLMSRSGCKVGDSSGTAAAHKVLPRGLSPARQLLLPASESPHWSGAPVKPSPQAAAVEVEEEDGSESEESAGPLLKGK PRALGGAAAG

靶点信息

研究背景	This gene encodes a member of the steroid receptor superfamily. The encoded protein mediates the physiological effects of progesterone, which plays a central role in reproductive events associated with the establishment and maintenance of pregnancy. This gene uses two distinct promoters and translation start sites in the first exon to produce several transcript variants, both protein coding and non-protein coding. Two of the isoforms (A and B) are identical except for an additional 165 amino acids found in the N-terminus of isoform B and mediate their own response genes and physiologic effects with little overlap.
基因ID	5241
基因名	PGR
Swiss	P06401
别名	NR3C3;PR;PGR

产品验证



Western blot analysis of Progesterone Receptor expressed in BT-474, Mouse brain, Rat brain using Progesterone Receptor Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn