

FBXO21 Rabbit pAb

货号: B13015

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	71kDa/72kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,MCF7,Mouse liver,Mouse testis,Rat brain
细胞定位	cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 440-621 of human FBX O21 (NP_055817.1).
序列	PEKVL DILQHIQTLDPGQHGA VGYLVQHTLEHIERKKEEVGVEVKLR SDEKHRDVCYSIGLIMKHKRYGYN CVIYGWDPTC MMGHEWIRNMNVHSLPHGHHQPFY NVLVEDGSCRYAAQENLEYNVEPQEISHPDVG RYFSEFTGTHYIPNAELEIRYPEDL EFVYETVQNIYSAKKENIDE

靶点信息

研究背景	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class. Alternative splicing of this gene generates 2 transcript variants.
基因ID	23014
基因名	FBXO21
Swiss	O94952
别名	FBXO21;FBX21

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn