

LASP1 Rabbit pAb

货号: B13038

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC IP
推荐浓度	WB: 1:1000 - 1:3000 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	23kDa/29kDa/36kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,BT-474,SKOV3,U-87MG,LO2,Mouse liver,Mouse brain
细胞定位	Cytoplasm,cell cortex,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-205 of human LASP 1 (NP_006139.1).
序列	RMGPSGGEGMEPERRDSQDGSSYRRPLEQQQPHHIPTSAPVYQQPQQQPVAQSYGGYKEPAAPVSIQRSAPGGGGK

靶点信息

研究背景	This gene encodes a member of a subfamily of LIM proteins, characterized by a LIM motif and a domain of Src homology region 3, and also a member of the nebulin family of actin-binding proteins. The encoded protein is a cAMP and cGMP dependent signaling protein and binds to the actin cytoskeleton at extensions of the cell membrane. The encoded protein has been linked to metastatic breast cancer, hematopoietic tumors such as B-cell lymphomas, and colorectal cancer.
------	---

基因ID	3927
基因名	LASP1
Swiss	Q14847
别名	LASP1;Lasp-1;MLN50

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn