

# TSR1 Rabbit pAb

货号: B13041

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>IHC:</b> Homo sapiens                                                                             |
| 应用    | <a href="#">WB</a>                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 91kDa                                                                                                |
| 实测分子量 | 110kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HepG2,Mouse testis                                                                                   |
| 细胞定位  | Nucleus,nucleolus                                                                                    |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                                                                               |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 455-804 of human TSR1 (NP_060598.3).                                                                                                                                                                                                                                                            |
| 序列   | HDDL YDKK VDEEAEAKMLEKYQERLEEMFPDEVDTPRDVAARIRFQKYRGLKSFRTSPWDPKENLPQDYARIFQFQNFT<br>NTRKSIFKEVEEKEVEGAEVGWYVTLHVSEVPVSVVECFRQGTPLIAFSLLPHEQKMSVLNMVRRDPGNTEPVKAKEELIF<br>HCGFRRFRASPLFSQHTAADKHKLQRFLTADMALVATVYAPITFPASVLLFKQKSNGMHSLIATGHLSVDPDRMVIKRV<br>VLSGHPFKIFTKMAVVRYMFFNREDVLWFKPVELRTKWGRRGHIKEPLGTHGHMKCSFDGKLKSQD TVLMNLYKRVFPK<br>WTYDPYVPEPVPWLKSEISSTVPQGGME |

靶点信息

|      |                                                                           |
|------|---------------------------------------------------------------------------|
| 研究背景 | Required during maturation of the 40S ribosomal subunit in the nucleolus. |
| 基因ID | 55720                                                                     |

|       |        |
|-------|--------|
| 基因名   | TSR1   |
| Swiss | Q2NL82 |
| 别名    | TSR1   |

产品验证



Western blot analysis of TSR1 expressed in HepG2,Mouse testis using TSR1 Rabbit pAb at 1:1000. Second ary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)