

WNT2 Rabbit pAb

货号: B13042

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	40kDa
实测分子量	43KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,NIH/3T3,Rat brain
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 26-360 of human WNT2 (NP_003382.1).
序列	SWWYMRATGGSSRVMCDNVPGLVSSQRQLCHRHPDVMRAISQGVAEWTAECQHQFRQHRWNCNTLDRDHSLFGRV LLRSSRESAFVYAISSAGVVFAITRACSQGEVKSCSDPKKMGS AKDSKGIFDWGGCSDNIDYGIKFARAFVDAKERKGK DARALMNLHNNRAGRKAVKRFLKQECKCHGVSGSCTLRTCWLAMADFRKTGDYLRKYNGAIQVVMNQDGTGFTVAN ERFKKPTKNDLVYFENSPDYCIRDREAGSLGTAGRVCNLT SRGMDSCEVMCCGRGYDTSHVTRMTKCGCKFWCCAVR CQDCLEALDVHTCKAPKNADWTTAT

靶点信息

研究背景	This gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Alternatively spliced transcript variants have been identified for this gene.
基因ID	7472
基因名	WNT2
Swiss	P09544
别名	WNT2;INT1L1;IRP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn