

HAS3 Rabbit pAb

货号: B13048

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	31kDa/62kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HL-60,SKOV3,SW480,Mouse heart,Rat liver,Rat heart
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 67-281 of human HAS3 (NP_619515.1).
序列	HRRMRRAGQALKLPSPRRGSVLALCIAAYQEDPDYLRKCLRSAQRISFPDLKVVMMVVDGNRQEDAYMLDIFHEVLGGTEQA GFFVWRSNFHEAGEGETEASLQEGMDRVRDVVRASTFSCIMQKWGGKREVMTAFKALGDSVDYIQVCDSDTVLDPA CTIEMLRVLEEDPQVGGVGGDVQPPGKGMAVEDDQVQAAQVRATEAWSVHQRHVSREQ

靶点信息

研究背景	The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyalu ronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, t his protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multi ple transcript variants.
------	--

基因ID	3038
基因名	HAS3
Swiss	O00219
别名	HAS3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn