

# FABP6 Rabbit pAb

货号: B13054

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                            |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>IF:</b> Homo sapiens                                                                                |
| 应用    | <a href="#">WB</a> <a href="#">IF/ICC</a>                                                              |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IF/ICC:</b> 1:50 - 1:200                                               |
| 理论分子量 | 14kDa/19kDa                                                                                            |
| 实测分子量 | 18kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | MCF7,A375                                                                                              |
| 细胞定位  | Cytoplasm                                                                                              |
| 纯化    | Affinity purification                                                                                  |

## 抗原信息


|      |                                                                                                                                     |
|------|-------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human FABP6 (NP_001436.1).                   |
| 序列   | MAFTGKFEMESEKNYDEFMKLLGISSDVIEKARNFKIVTEVQQDGQDFTWSQHYSGGHTMTNKFTVGKESNIQTMGGKT<br>FKATVQMEGGKLVNFPNYHQTSEIVGDKLVEVSTIGGVTYERVSKRLA |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes the ileal fatty acid binding protein. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. FABP6 and FABP1 (the liver fatty acid binding protein) are also able to bind bile acids. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism. Transcript variants generated by alternate transcription promoters and/or alternate splicing have been found for this gene. |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                                  |
|-------|--------------------------------------------------|
| 基因ID  | 2172                                             |
| 基因名   | FABP6                                            |
| Swiss | P51161                                           |
| 别名    | FABP6;I-15P;I-BABP;I-BALB;I-BAP;ILBP;ILBP3;ILLBP |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)