

TUSC2 Rabbit pAb

货号: B13057

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Sea cucumber
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	12kDa
实测分子量	12kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human TUSC2 (NP_009206.1).
序列	MGASGSKARGLWPFASAAGGGGSEAAGAEQALVRPRGRAVPPFVFTRRGSMFYDEDDGLAHEFYEETIVTKNGQKRAK LRRVHKNLIPQGIVKLDHPRIHVDFPVILYEV

靶点信息

研究背景	This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.
基因ID	11334

基因名	TUSC2
Swiss	O75896
别名	TUSC2;C3orf11;FUS1;PAP;PDAP2

产品验证



Western blot analysis of TUSC2 expressed in human esophagus cancer,mouse liver,rat lung using TUSC2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn