

EGR3 Rabbit pAb

货号: B13066

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	38kDa/42kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EGR3 (NP_004421.2).
序列	MTGKLAEKLPVTMSSLLNQLPDNLYPEEIPSA N LFSGSSDSVVHYNQMATENVMDIGLTNEKPNPELSYSGSFQPAPGNK TVTYLGKFAFDSPSNWCQDNIISLMSAGILGVPPASGALSTQTSTASMVQPPQGDVEAMYPALPPYSNCGDLYSEPVSFH DPQGNPGLAYSPQDYQSAK P ALDSNLFPMIPDYNLYHHPNDMGSIPEHKPFQGM D PIRVNPPPITPLETIKAFKDKQIH P G FGSLPQPP

靶点信息

研究背景	This gene encodes a transcriptional regulator that belongs to the EGR family of C2H2-type zinc-finger proteins. It is an immediate-early growth response gene which is induced by mitogenic stimulation. The protein encoded by this gene participates in the transcriptional regulation of genes in controlling biological rhythm. It may also play a role in a wide variety of processes including muscle development, lymphocyte development, endothelial cell growth and migration, and neuronal development. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
基因ID	1960
基因名	EGR3
Swiss	Q06889
别名	EGR3;EGR-3;PILOT

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn