

# MyD88 Rabbit pAb

货号: B13082

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Rat                                                                                            |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>IHC:</b> Rattus norvegicus                                                                        |
| 应用    | <a href="#">WB</a> <a href="#">IHC</a>                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IHC:</b> 1:50 - 1:200                                                |
| 理论分子量 | 15kDa/20kDa/28kDa/31kDa/33kDa/34kDa                                                                  |
| 实测分子量 | 33KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | A-549,Rat lung                                                                                       |
| 细胞定位  | Cytoplasm                                                                                            |
| 纯化    | Affinity purification                                                                                |


抗原信息

|      |                                                                                                                                                                                                                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-296 of human MyD88 (Q99836).                                                                                                                                                                                                       |
| 序列   | MAAGGPGAGSAAPVSSTSSLPLAALNMRVRRRLSLFLNVRTQVAADWTALAEEMDFEYLEIRQLETQADPTGRLLDAWQ<br>GRPGASVGRLLELLTKLGRDDVLLELGPSIEEDCQKYILKQQQEAEKPLQVAAVDSSVPRTAELAGITTLDDPLGHMPERF<br>DAFICYCPSDIQFVQEMIRQLEQTNYRLKLCVSDRDVLPGTCVWSIASELIEKRCRRMVVVVSDDYLSKECDFQTKFALSL<br>SPGAHQKRLIPIKYKAMKKEFPSILRFITVCDYTNPCTKSWFWTRLAKALSLP |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene encodes a cytosolic adapter protein that plays a central role in the innate and adaptive immune response. This protein functions as an essential signal transducer in the interleukin-1 and Toll-like receptor signaling pathways. These pathways regulate that activation of numerous proinflammatory genes. The encoded protein consists of an N-terminal death domain and a C-terminal Toll-interleukin1 receptor domain. Patients with defects in this gene have an increased susceptibility to pyogenic bacterial infections. Alternate splicing results in multiple transcript variants. |
| 基因ID  | 4615                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 基因名   | MYD88                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Swiss | Q99836                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 别名    | MYD88;MYD88D;myeloid differentiation primary response 88;MyD88                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)