

LTBP3 Rabbit pAb

货号: **B13093**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	134kDa/139kDa
实测分子量	145kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1004-1303 of human LTBP3 (NP_001123616.1).
序列	ICKEGKCVNTQPGYECYCKQGFYYDGNLLECVDVDECLDESNCRNGVCENTRGGYRCACPPAEYSPAQRQCLSPEEMDVDECQDPAACRPGRCVNLPGSYRCECRPPWVPGPSGRDCQLPESPAERAPERRDVCWSQRGEDGMCAGPLAGPALTFDDCCCRQGRGWGAQCRPCPPRGAGSHCPTSQSESNSFWDTSPLLL GKPPRDEDSSEEDSDECRCVSGRCVPRPGGAVCECPGGFQLDASRARCVDIDECRELNQRGLLCKSERCVNTSGSFRCVCKAGFARSRPHGACVPQRRR

靶点信息

研究背景	The protein encoded by this gene forms a complex with transforming growth factor beta (TGF-beta) proteins and may be involved in their subcellular localization. Activation of this complex requires removal of the encoded binding protein. This protein also may play a structural role in the extracellular matrix. Three transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	4054
基因名	LTBP3
Swiss	Q9NS15
别名	LTBP3;DASS;LTBP-3;LTBP2;STHAG6;pp6425

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn