

RNF149 Rabbit pAb

货号: B13096

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	43kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,mouse pancreas,rat brain
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 221-400 of human RNF149 (NP_775918.2).
序列	FYYIQRFLYTGSQIGSQSHRKETKKVIGQLLLHTVKHGEKGIDVDAENCAVCIENFKVKDIIRILPCKHIFHRICIDPWLLDHRTCPMCKLDVIKALGYWGEPGDVQEMPAPESPGRDPAANLSLALPDDDGSDSSPPSASPASEPQCDPSFKGDAGENTALLEAGRSDSRHGGPIS

靶点信息

研究背景	E3 ubiquitin-protein ligase. Ubiquitinates BRAF, inducing its proteasomal degradation.
基因ID	284996

基因名	RNF149
Swiss	Q8NC42
别名	RNF149;DNAPTP2

产品验证



Western blot analysis of RNF149 expressed in mouse pancreas, rat brain, HeLa using RNF149 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn