

# HCK Rabbit pAb

货号: B13109

产品信息

|       |                                                                                                                                                                                                                            |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                                                                                                                            |
| 宿主    | Rabbit                                                                                                                                                                                                                     |
| 克隆性   | Polyclonal                                                                                                                                                                                                                 |
| 预测反应  | <b>WB:</b> Mus musculus                                                                                                                                                                                                    |
| 应用    | <a href="#">WB</a>                                                                                                                                                                                                         |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                                                                                                                                                  |
| 理论分子量 | 57kDa/59kDa                                                                                                                                                                                                                |
| 实测分子量 | 55-60kDa                                                                                                                                                                                                                   |
| 形式    | Liquid                                                                                                                                                                                                                     |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.                                                                                                                       |
| 偶联物   | Unconjugated                                                                                                                                                                                                               |
| 阳性对照  | Raji,Mouse spleen,Mouse lung,Rat spleen                                                                                                                                                                                    |
| 细胞定位  | Cell junction,Cell membrane,Cell projection,Cytoplasm,Cytoplasmic vesicle,Golgi apparatus,Lipid-anchor, Lysosome,Membrane,Nucleus,caveola,cytoskeleton,cytosol,cytosol,focal adhesion,podosome membrane, secretory vesicle |
| 纯化    | Affinity purification                                                                                                                                                                                                      |


抗原信息

|      |                                                                                                                                                                                               |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 324-505 of human HCK (NP_001165600.1).                                                                          |
| 序列   | LKSDGSKQPLPKLIDFSAQIAEGMAFIEQRNYIHRDLRAANILVSASLVCKIADFGLARVIEDNEYTAREGAKFPIKWTAPEA<br>INFGSFTIKSDVWSFGILLMEIVTYGRIPYPGMSNPEVIRALERGYRMPRPENCPEELYNIMMRCWKNRPEERPTFEYIQSVL<br>DDFYTATESQYQQQP |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | The protein encoded by this gene is a member of the Src family of tyrosine kinases. This protein is primarily hemopoietic, particularly in cells of the myeloid and B-lymphoid lineages. It may help couple the Fc receptor to the activation of the respiratory burst. In addition, it may play a role in neutrophil migration and in the degranulation of neutrophils. Multiple isoforms with different subcellular distributions are produced due to both alternative splicing and the use of alternative translation initiation codons, including a non-AUG (CUG) codon. |
| 基因ID  | 3055                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 基因名   | HCK                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Swiss | P08631                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 别名    | HCK;JTK9;p59Hck;p61Hck                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)