

ZNF582 Rabbit pAb

货号: B13110

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	60kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-180 of human ZNF582 (NP_653291.1).
序列	SNLVSGLAVSKPDVISFLEQGKEPWMVERVVSGGLCPVLESRYDTKELFPKQHVVYEVESPQWEIMESLTSYGLECSSFQDDWECRNQFDRQQGNPDRHFHQMIIRHEEMPTFDQHASLTFYQKIHTREKPFGYNKCRKDF

靶点信息

研究背景	The protein encoded by this gene is a zing finger protein and putative transcription factor that is highly methylated in cervical cancers. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	147948

基因名	ZNF582
Swiss	Q96NG8
别名	ZNF582

产品验证



Western blot analysis of ZNF582 expressed in Mouse brain using ZNF582 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn