

RIN1 Rabbit pAb

货号: B13119

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	77kDa/84kDa
实测分子量	84kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse testis,Rat brain
细胞定位	Cytoplasm,Membrane,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RIN1 (N P_004283.2).
序列	MESPGESGAGSPGAPSPSSFTTGHLAREKPAQDPLYDVPNASGGQAGGPQRPGRVVSLRERLLLTRPVWLQLQANAAA ALHMLRTEPPGTFLVRKSNTRQCQALCMRLPEASGPSFVSSHYLESPGGVSLEGSELMFPDLVQLICAYCHTRDILLPLQL PRAIHHAATHKELEAISHLGIEFWSSSLNIKAQRGPAGGPVLPQLKARSPQELDQGTGAALCFFNPLFPGDLGPTKREKFKR SFKVRVSTETSSPLSPPAVPPPPVPVLPGAVPSQTER

靶点信息

研究背景	Ras effector protein, which may serve as an inhibitory modulator of neuronal plasticity in aversive memor y formation. Can affect Ras signaling at different levels. First, by competing with RAF1 protein for binding to activated Ras. Second, by enhancing signaling from ABL1 and ABL2, which regulate cytoskeletal remod eling. Third, by activating RAB5A, possibly by functioning as a guanine nucleotide exchange factor (GEF f or RAB5A, by exchanging bound GDP for free GTP, and facilitating Ras-activated receptor endocytosis.
------	--

基因ID	9610
基因名	RIN1
Swiss	Q13671
别名	RIN1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn